

ADC Photoelectric Timing Adjustment

1. Adjustment tool.....	2
2. Preparation	2
3. Photoelectric ADC Timing Adjustment	2

1. Adjustment tool

N/A

2. Preparation

N/A

* Since TM is used, maintenance key is not used (remain to "OFF").



3. Photoelectric ADC Timing Adjustment

(1) Select [Maintenance] > [Service] > [ADC Photoelectric Timing Adjustment] > [Select].

The screenshot displays the c703 service menu. At the top, it shows 'Stand By' and 'User ID HHT'. The main menu has tabs for 'Settings', 'System', 'Application', 'Special Wash', and 'Setup'. The 'Maintenance' tab is active, showing a list of tasks. The 'Service' option is highlighted in red. The 'ADC Photoelectric Timing Adjustment' task is also highlighted in red, and the 'Select' button next to it is also highlighted in red. The interface includes a right sidebar with navigation icons for 'Pre-Routine', 'Overview', 'Maintenance', 'Alarms', and 'Report'. The bottom status bar shows various system indicators like 'Routine', 'Sample Status', 'Reagents', 'Calibration', 'QC', and 'Data Monitor'.

(2) Select [Perform].

The screenshot shows a software interface for maintenance. The main window is titled 'Maintenance' and contains a sub-window titled 'ADC Photoelectric Timing Adjustment'. Inside this sub-window, there are three icons representing different components: 'SSU', 'SB1', and 'c7-1'. A 'Details' button is located in the top right corner of the sub-window. At the bottom of the sub-window, there is a text box labeled 'Comment on Report'. Below the sub-window, there are two buttons: 'Perform' and 'Cancel'. The 'Perform' button is highlighted with a red rectangular box.

(3) Confirm that no alarm occurs during TM operation (takes about 30 minutes).

If alarm occurs, check the LED lamp, wiring/connector (LED lamp - CDIO8, photometer - CLADC - NCOMCAD200) and Reaction cell. Replace if necessary.

(Refer to "S07-01 Alarm Code" for details on alarms.)

End of Document.